



La Radiologie Vétérinaire,  
Notre Métier.

# Option : Swivelling monobloc



The 'monobloc' is the part of the system at the top of the column and that contains the Rx tube.

The primary X ray field can be swivelled upto 45° to the right or to the left and enables to take radiographs of specific incidences such sinusitis for example.

The 'zero' or rest position is marked by a notch.

A handle is added to the monobloc to ease the swivelling.

**GeR-International**

1 rue de Bourgogne 89140 CUY France

Standard : 33 (0) 3-86-64-07-31 Fax : 33 (0) 3-86-65-90-03

E-mail : [contact@ger-int.com](mailto:contact@ger-int.com) Site Web : [www.ger-int.com](http://www.ger-int.com)