



La radiologie vétérinaire notre métier

CONTENTION ACCESSORIES



This options enables to tightly maintain and position the animal on the table top in order to perform quality images.

The option is composed of 2 elements :

- The table top features a specific aluminium profile along the panel.
- The contetion system itself.

Characteristics :

- Le contention system moves all the way along the table panel.
- The fabric (165mm width) is stretched thanks to a rolling knob, is easily removeable and is washable with laundry machines.
- A simple press on a button releases the tension of the system.
- The contention systems can be completely removed from the table top.

GeR-International : 1 rue de Bourgogne 89140 CUY France
Standard : 33 (0) 3-86-64-07-31 Fax : 33 (0) 3-86-65-90-03
E-mail : ger2@wanadoo.fr Site Web : www.ger-int.com